

MODEL SFPvh 6501 PERFECTION

SOLID DOOR LAB FREEZER -9°C TO -35°C

Capacity [Cu. Ft.]	29.1 Cubic Feet
Capacity [Liters]	831 Liters
Temperature Range [°C]	-10°C to -35°C
Temperature Range [°F]	+14°F to -31°F

Dimensions	Width x Depth x Height
Exterior [Inches]	27.5" x 34.25" x 84.25"
Interior [Inches]	23.75" x 26" x 59.25"
Exterior [Centimeters]	70 x 87 x 214
Interior [Centimeters]	60 x 66 x 150

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R290
Exterior	Steel, Painted White
Interior	Stainless Steel
Insulation	Polyurethane Foam
Interior Lighting	LED Ceiling Light
Self-Closing Door	1 Door, Automatic, Recessed
Lockable Door	Keyed Lock
Door Material	Insulated Steel
Porthole	1" Diameter
Casters	2 Swivel Locking, 2 Fixed
Shelving	4 Adjustable Shelves
Shelving Weight Limit	132 Lbs.
Approximate KJ/Hour	875
Approximate Decibels	60 dB

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 10A
NEMA Configuration	NEMA 5-15
Certification(s)	CSA Listed Energy Star Rated

Control:

Display Type	Microprocessor, Digital
Alarms	High / Low Temp
Min / Max Temperature	Display and Reset
Alarm Signal	Audible and Visual
Alarm Relay	Dry Contacts
Alarm Forwarding	Yes
Data Logger	Integrated, USB Interface



Warranty 3 Years Parts & Labor