

MODEL U85-11

UPRIGHT ULTRA LOW FREEZER 0°C TO -40°C

Capacity [Cu. Ft.]	11 Cu. Ft.
Capacity [Liters]	311 Liters
Temperature Range [°C]	-40°C to -85°C
Temperature Range [°F]	-40°F to -120°F

Dimensions Width x Depth x Height

Exterior [Inches] 28" x 34.5" x 79.5"

Interior [Inches] 17" x 22" x 51.25"

Exterior [Centimeters] 71 x 87 x 201

Interior [Centimeters] 43 x 55 x 130

Refrigeration System Cascade-Style System

Compressors 2, Hermetic Compressors

Air Circulation Gravity Flow

Defrost Method Manual

Refrigerant HCFC & CFC Free

Exterior Powder Coated Steel

Interior Zinc Coated Galvanized Steel

Insulation Polyurethane (Ecomate) Foam

Inner Doors 5, Stainless Steel

Lockable Door Lockable with Padlock

Porthole Optional 1/2"

Casters 4, Swivel with Locks

Storage 4 Adjustable Shelves, with

1 Base Shelf

Shipping Weight 850 Lbs.

Electrical

Standard Voltage 115 Volts, 60 Hertz, 1 Phase

Amps / Run Amps / Plug 20A | 7.3A | NEMA 5-20

(Optional Voltages) 208 Volts, 230 Volts

Certification(s) ETL Listed

Controller

Display Type Microprocessor, Digital

Alarms High/Low Temp, Power

Failure, Low Battery

Alarm Relay Yes, Dry Contacts

Automatic Restart Yes, Automatically Restarts

After Power Failure

Backup Battery Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor