


ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW FREEZER RACKS

FOR ALL STANDARD 2.00 & 3.00-INCH STORAGE BOXES

For All Models of Chest Freezers and Refrigerators

The highest quality racks, constructed from corrosion resistant stainless steel, including handles. Racks listed are for standard 2.00 or 3.00-inch high cardboard, aluminum stainless steel, or plastic boxes. Custom racks available upon request.



CHEST RACK CAPACITY:

- 3 cu. ft. chest can hold 8 racks
- 5 cu. ft. chest can hold 15 racks
- 9 cu. ft. chest can hold 16 racks
- 12 cu. ft. chest can hold 27 racks
- 14 cu. ft. chest can hold 30 racks
- 17 cu. ft. chest can hold 36 racks
- 22 cu. ft. chest can hold 36 racks
- 27 cu. ft. chest can hold 48 racks

UPRIGHT RACK CAPACITY:

- 13 cu. ft. upright can hold 15 racks
- 18 cu. ft. upright can hold 20 racks
- 22 cu. ft. upright can hold 25 racks
- 25 cu. ft. upright can hold 30 racks
- 30 cu. ft. upright can hold 35 racks

! All racks come standard with boxes and 81 cell dividers.

2" Box Dimensions (Inches)	
W	5.25"
D	5.25"
H	2"

3" Box Dimensions (Inches)	
W	5.25"
D	5.25"
H	3"

Model Number	Type	Length	Height	Depth	Box Size	Box Capacity
16-5-3	Chest	5.62"	15.87"	5.5"	3" Boxes	5 Boxes
16-7-2	Chest	5.62"	15.43"	5.5"	2" Boxes	7 Boxes
19-6-3	Chest	5.62"	19.12"	5.5"	3" Boxes	6 Boxes
19-9-2	Chest	5.62"	19.81"	5.5"	2" Boxes	9 Boxes
29-9-3	Chest	5.62"	29.25"	5.5"	3" Boxes	9 Boxes
29-13-2	Chest	5.62"	28.25"	5.5"	2" Boxes	13 Boxes
17-9-3	Upright	16.31"	9.37"	5.5"	3" Boxes	9 Boxes
17-12-2	Upright	16.31"	8.87"	5.5"	2" Boxes	12 Boxes
17-12-2-DR	Upright (Drawer)	16.5"	9.43"	5.5"	2" Boxes	12 Boxes
23-12-3	Upright	21.62"	9.37"	5.5"	3" Boxes	12 Boxes
23-16-2	Upright	21.62"	8.87"	5.5"	2" Boxes	16 Boxes
23-16-2-DR	Upright (Drawer)	22"	9.43"	5.5"	2" Boxes	16 Boxes