

SPECIAL PRICE – SCRATCH & DENT

SO-LOW CHEST STYLE FREEZER TO -40°C MODEL CH40-9-SA

9 Cubic Ft. (255 liters) Chest – Temperature Range: 0°C to -40°C

Specifications:

*Capacity	9 cubic ft. / 255 liters
*Temperature Range	0°C to -40°C
*Dimensions	
- Exterior w x d x h	39" x 23.5" x 33.5"
- Interior w x d x h	33.5" x 17.5" x 27.50"

Special Price: \$2,995.00
Model CH40-9

CONTROL

Digital Temperature control with two Temperature displays. One display shows the set temperature, the other shows the actual temperature in the freezer.

ALARM

Battery operated alarm system provides An audible and visual signal in a high/low temperature alarm condition. The alarm will also signal in the event of a power failure. Also provided is an alarm silencing and battery test feature.

MANUAL DEFROST

CASTERS FOR EASY MOBILITY

1/3HP HERMETIC MOTOR COMPRESSOR

115 VOLT, 60 HERTZ, 1 PHASE

WARRANTY One Year Parts
90 Days Labor

TECHNICAL SUPPORT

Supplied by Nation-Wide-Network of service companies.

ETL Electrical Approval



SPECIAL OFFER # SA25