

SPECIAL PRICE – OVERSTOCK

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-10GD-SA

10 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

* Capacity 10 cubic ft.
* Temperature Range 2°C to 8°C / 36°F to 46°F
* Dimensions

* Exterior w x d x h 24" x 26" x 65.5"
60 x 66 x 166 (cm)

* Interior w x d x h 20.75" x 21" x 37.75"
52 x 53 x 96 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Door	Single glass, locking
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four locking
Shelves	Three Adjustable
Interior Lighting	LED Light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Compressor	1/5 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	2.0 / 15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	300 lbs.
Warranty-(Domestic)	One year parts and labor.

Special Price: \$2,910.00



SPECIAL OFFER # SA26