

SPECIAL PRICE – SCRATCH & DENT

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-37GD-SA

30 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

- * Capacity 37 cubic ft. / 1047 liters
- * Temperature Range 2°C to 8°C / 36°F to 46°F

Special Price: \$4,370.00

When ordering please specify Celsius or Fahrenheit Setting

- * Dimensions
- * Exterior w x d x h 47.12" x 30.12" x 79"
- * Interior w x d x h 43 x 22.25" x 60.5"

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter
Interior Finish	White coated steel
Exterior Finish	White coated steel
Doors	Double glass – Lockable
Casters	Four, locking
Shelves	Ten, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane
Refrigerant:	R290A
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps	8.3/15
Supply Plug	5-15 NEMA UL listed
Shipping Weight	500 lbs.
Warranty-(Domestic)	One year parts and labor.



SPECIAL OFFER # SA27