

SPECIAL PRICE – SCRATCH & DENT

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-30SGD-SA

30 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

Special Price:

\$4,490.00

- * Capacity 30 cubic ft. / 849 liters
- * Temperature Range 2°C to 8°C / 36°F to 46°F

When ordering please specify Celsius or Fahrenheit Setting

- * Dimensions
- * Exterior w x d x h 47.25" x 30.25" x 81"
120 x 76 x 205 (cm)
- * Interior w x d x h 43 x 22" x 54"
109 x 55 x 137 (cm)

Features:

- | | |
|---------------------|---|
| Temperature Control | Digital display
With Hi/Low Alarm
Alarm Relay Dry Contacts
Min/Max Memory
Power Failure Relay |
| Access Port | 1/2" Diameter |
| Interior Finish | White coated steel |
| Exterior Finish | White coated steel |
| Doors | Double glass – sliding
Lockable |
| Casters | Four, locking |
| Shelves | Ten, Adjustable |
| Interior Lighting | LED light |
| Insulation | Polyurethane |
| Refrigerant: | R290A |
| Air Circulation | Forced air |
| Defrost | Cycle Defrost |
| Voltage: Domestic | 115/60/1 |
| Running Amps | 8.3/15 |
| Supply Plug | 5-15 NEMA
UL listed |
| Shipping Weight | 500 lbs. |
| Warranty-(Domestic) | One year parts and labor. |



SPECIAL OFFER # SA28