

SPECIAL PRICE – SCRATCH & DENT

SO-LOW UNDERCOUNTER FREEZER MODEL MV40-4UCF-SA

4 Cubic Ft. Upright – Temperature Range: -30°C to -40°C

Specifications:

* Capacity	4 cubic ft.
* Temperature Range	-30°C to -40°C
* Dimensions	
* Exterior w x d x h	24" x 24" x 33" 61 x 61 x 83 (cm)
* Interior w x d x h	17.5" x 15.5" x 24.75" 44 x 39 x 62 (cm)

Special Price: \$2,450.00

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts
Access Port	3/8" side
Door	Single solid, locking
Shelves	Three adjustable
Insulation	Polyurethane foam
Refrigerant:	R290
Compressor	1/4 hp
Defrost	Manual
Voltage: Domestic	115/60/1
Running Amps/Breaker	5/15
Supply Plug	5-15P NEMA
Shipping Weight	150 lbs.
Warranty-(Domestic)	Two year parts and labor



SPECIAL OFFER # SA31