

SPECIAL PRICE – SCRATCH & DENT

SO-LOW LAB & PHARMACY FREEZER MODEL DHS25-25SD

25 Cubic Ft. (679 liters) Solid Door Upright

Temperature Range: -18°C to -25°C / 0°F to -13°F

Specifications:

Special Price: \$6,520.00

| | |
|----------------------|-------------------------------|
| * Capacity | 25 cubic ft. / 679 liters |
| * Temperature Range | -18°C to -25°C / 0°F to -13°F |
| * Dimensions | |
| * Exterior w x d x h | 28.35" x 31.5" x 79.25" |
| * Interior w x d x h | 23.62" x 26.18" x 59.84" |

Features:

| | |
|----------------------|--|
| Temperature Control | Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Power Failure Relay Battery back up |
| Access Port | 1/2" Diameter |
| Door | Single solid – lockable |
| Casters | heavy duty, swivel |
| Shelves | Three Sliding Adjustable |
| Interior Lights | Led light |
| Insulation | Polyurethane CFC Free |
| Refrigerant: | R-290 |
| Air Circulation | Forced air |
| Evaporator | Fin and tube, forced air |
| Thermostat | Adjustable |
| Defrost | Automatic |
| Voltage: | 115/60/1 |
| Running Amps/Breaker | 10.3/15 |
| Supply Plug | 5-15 NEMA |
| Shipping Weight | 355 lbs. |
| Warranty-(Domestic) | One Year Parts. |



Specifications subject to change without notice.

SPECIAL OFFER # SA34