

SPECIAL PRICE – SCRATCH & DENT

SO-LOW LAB & PHARMACY FREEZER

MODEL DHS25-25SD

25 Cubic Ft. (679 liters) Solid Door Upright

Temperature Range: -18°C to -25°C / 0°F to -13°F

Specifications:

* Capacity	25 cubic ft. / 679 liters
* Temperature Range	-18°C to -25°C / 0°F to -13°F
* Dimensions	
* Exterior w x d x h	28.35" x 31.5" x 79.25"
* Interior w x d x h	23.62" x 26.18" x 59.84"

Special Price: \$5,500.00

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Power Failure Relay Battery back up
Access Port	1/2" Diameter
Door	Single solid – lockable
Casters	heavy duty, swivel
Shelves	Three Sliding Adjustable
Interior Lights	Led light
Insulation	Polyurethane CFC Free
Refrigerant:	R-290
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Automatic
Voltage:	115/60/1
Running Amps/Breaker	10.3/15
Supply Plug	5-15 NEMA
Shipping Weight	355 lbs.
Warranty-(Domestic)	One Year Parts.



Specifications subject to change without notice.

SPECIAL OFFER # SA34