

## SPECIAL PRICE – SCRATCH & DENT

### SO-LOW LAB & PHARMACY REFRIGERATOR

#### MODEL DHS4-30SD

**30 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F**

**Solid Door**

#### Specifications:

* Capacity	30 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	33.85" x 37.80" x 81.12"
* Interior w x d x h	27.50" x 30.80" x 61.84"

#### Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Power Failure Relay Battery back up
Door	Single solid, locking
Access Port	1/2" Diameter
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four locking
Shelves	Four Sliding Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	6.0 / 15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	420 lbs.
Warranty-(Domestic)	One Year Parts.

**Special Price: \$5,340.00**



Specifications subject to change without notice

## SPECIAL OFFER # SA36