

## SPECIAL PRICE – SCRATCH & DENT

### SO-LOW LAB & PHARMACY FREEZER MODEL DHS25-49SD

49 Cubic Ft. (1387 liters) Upright - Temperature Range: -18°C to -25°C / 0°F to -13°F  
2 Solid Doors

#### Specifications:

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	-18°C to -25°C / 0°F to -13°F
* Dimensions	
* Exterior w x d x h	56.69" x 31.5" x 79.25"
* Interior w x d x h	51.97" x 26.18" x 59.84"

#### Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Power Failure Relay Battery back up
Access Port	1/2" Diameter
Doors	Two doors solid - lockable
Casters	Heavy duty, swivel
Shelves	Six Sliding Adjustable
Interior Lights	Led lights
Insulation	Polyurethane CFC Free
Refrigerant:	R-290
Compressor	1 hp
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Automatic
Voltage:	115/60/1
Running Amps/Breaker	6.2/20
Supply Plug	5-20 NEMA
Shipping Weight	450 lbs.
Warranty-(Domestic)	One year parts.

Special Price: **\$8,600.00**



Specifications subject to change without notice.

## SPECIAL OFFER # SA39