

SPECIAL PRICE – SCRATCH & DENT

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-8GD

8 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Glass Door

Specifications:

Special Price: \$2,400.00

* Capacity 8 cubic ft.
* Temperature Range 2°C to 8°C / 36°F to 46°F
When ordering please specify Celsius or Fahrenheit Setting

* Dimensions
* Exterior w x d x h 19" x 21" x 74.5"
48 x 53 x 189 (cm)
* Interior w x d x h 15.87" x 14.5" x 53"
40 x 36 x 134 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Door	Single glass, locking
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four locking
Shelves	Three Adjustable
Interior Lighting	LED Light
Insulation	Polyurethane foam
Refrigerant:	R600 CFC free
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	2.0 / 15
Supply Plug	5-15P NEMA
Shipping Weight	300 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice

SPECIAL OFFER # SA44