

## MODEL DH4-8GD | SA44

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

**Capacity [ Cu. Ft. ]** 6.17 Cu. Ft.  
**Capacity [ Liters ]** 175 Liters  
**Temperature Range [ °C ]** 2°C to 8°C  
**Temperature Range [ °F ]** 36°F to 46°F

**Dimensions** Width x Depth x Height  
**Exterior [ Inches ]** 19" x 20.75" x 71.5"  
**Interior [ Inches ]** 15.75" x 14.5" x 46.75"  
**Exterior [ Centimeters ]** 48 x 52 x 181  
**Interior [ Centimeters ]** 40 x 36 x 118

**Compressor** Hermetic Compressor  
**Air Circulation** Forced Air  
**Defrost Method** Cycle  
**Refrigerant** R600a  
**Exterior** Painted Steel  
**Interior** Painted Steel  
**Insulation** Polyurethane Foam  
**Lockable Door** 1 Door, Keyed Lock  
**Door Material** Glass  
**Interior Light** Adjustable LED  
**Porthole** 1/2" Diameter ( Top Back Wall )  
**Casters** 4  
**Storage** 5 Adjustable Shelves  
**Shipping Weight** 300 lbs.

#### Electrical

**Voltage:** 115 Volts 60 Hertz 1 Phase  
**Amperage Line | Running** 15 Amp Dedicated | 2.0A  
**NEMA Configuration** NEMA 5-15  
**Certification(s)** UL Listed  
Energy Star Rated

#### Controller

**Display Type** Microprocessor, Digital  
**Alarms** High/Low Temp  
**Min / Max Temperature** Display and Reset  
**Alarm Relay** Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 1 Year Parts & Labor.