

SPECIAL PRICE – OVERSTOCK

SO-LOW LAB & PHARMACY FREEZER MODEL DH30-23SDSS

23 Cubic Ft. Upright – Temperature Range: -18°C to -25°C / 0°F to -13°F

Specifications:

Special Price: \$4,650.00

- * Capacity 23 cubic ft.
- * Temperature Range -18°C to -20°C / 0°F to -4° F
- *When ordering please specify Celsius or Fahrenheit Setting*

- * Dimensions
- * Exterior w x d x h 27.25" x 33.75" x 81"
69 x 85 x 206 (cm)
- * Interior w x d x h 24" x 28.5" x 55"
61 x 72 x 140 (cm)

Features:

- Temperature Control Digital display
W/ Hi-Low Alarm
Alarm Relay Dry Contacts
Min/Max Memory
Power Failure Relay
- Access Port 1/2" Diameter Top Back Wall
- Doors Single Solid, locking
- Interior Finish Stainless Steel
- Exterior Finish Stainless Steel
- Casters Four
- Shelves Three
- Interior Lighting One 20-watt fluorescent
- Insulation Urethane foam
- Refrigerant: R404A CFC free
- Compressor 1/2 hp
- Air Circulation Forced air
- Evaporator Fin and tube, forced air
- Defrost Automatic
- Voltage: Domestic 115/60/1
- Running Amps/Breaker 10.9/15
- Supply Plug 5-15P NEMA
UL listed
- Shipping Weight 440 lbs.
- Warranty-(Domestic) 90 days parts and labor. Compressor parts additional 2 years, no labor.



Specifications subject to change without notice

SPECIAL OFFER # SA46