

## SPECIAL PRICE – OVERSTOCK

### SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-50PT

50 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

#### Specifications:

Special Price: \$12,750.00

- \* Capacity 50 cubic ft
- \* Temperature Range 2°C to 8°C / 36°F to 46°F
- \* **When ordering please specify Celsius or Fahrenheit Setting\***

- \* Dimensions
- \* Exterior w x d x h 52" x 37.5" x 81.5"  
132 x 95 x 207 (cm)
- \* Interior w x d x h 48" x 30" x 59"  
121 x 76 x 149 (cm)

#### Features:

- Temperature Control Digital display with  
Hi/Low alarm  
Alarm Relay Dry Contacts  
Min/Max Memory  
Power Failure Relay
- Interior Finish Stainless Steel
- Exterior Finish Stainless Steel
- Door Glass doors front and back, lockable
- Casters Four 3" swivel, locking
- Shelves Six, Adjustable
- Interior Light One
- Insulation: Polyurethane foam
- Refrigerant: R-290A
- Compressor 3/4 hp
- Air Circulation Forced air
- Defrost Cycle Defrost
- Voltage: 115/60/1
- Running Amps 8.0
- Supply Plug 5-15P NEMA  
UL listed
- Shipping Weight 650 lbs.
- Warranty-(Domestic) 1 year parts and labor. Compressor parts additional 2 years, no labor.



Specifications subject to change without notice.

## SPECIAL OFFER # SA48