

SPECIAL PRICE – SCRATCH & DENT

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHS4-25SD-SS

25 Cubic Ft. Upright –a Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

Special Price: **\$4,600.00**

* Capacity	25 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	28.35" x 31.5" x 79.25"
* Interior w x d x h	23.62" x 26.18" x 59.84"

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Power Failure Relay Battery back up
Door	Single solid, locking
Access Port	1/2" Diameter Top Back Wall
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four locking
Shelves	Four Sliding Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	6.0 / 15
Supply Plug	5-15P NEMA
Shipping Weight	400 lbs.
Warranty-(Domestic)	One year parts and labor Compressors additional 3 years, no labor.



SPECIAL OFFER # SA50