

## SPECIAL PRICE – OVERSTOCK

### SO-LOW UNDERCOUNTER FREEZER MODEL MV85-2

2 Cubic Ft. (54 liters) Upright – Temperature Range: -40°C to -85°C

#### Specifications:

Special Price: **\$5,500.00**

* Capacity	2 cubic ft. / 54 liters
* Temperature Range	-40°C to -85°C
* Dimensions	
* Exterior w x d x h	23" x 27" x 32" 60 x 70 x 81 (cm)
* Interior w x d x h	14" x 18" x 12.5" 36 x 46 x 32 (cm)

#### Features:

Control	Digital Datalogging Battery Back-up Hi/Low Alarm
Clearance	4" on sides and back
Doors	Single Solid, locking
Legs	Four, leveling
Insulation	Polyurethane CFC Free
Evaporator	Cold wall
Defrost	Manual
<b>Voltage: Domestic</b>	<b>220-50/60-1 – 6.3 amps</b>
Supply Plug	NEMA 5-15P ETL listed
Shipping Weight	265 lbs.
Warranty	1 Year Parts. 90 Days Labor.



Specifications subject to change without notice

## SPECIAL OFFER # SA51