

## *SO-LOW LAB CHROMATOGRAPHY PHARMACY REFRIGERATOR MODEL DHS4-49GD-CH*

**49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F**

**Two Swinging Glass Doors**

### **Specifications:**

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	56.69" x 33.5" x 79.25"
* Interior w x d x h	51.97" x 26.18" x 59.84"

### **Features:**

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Power Failure Relay Battery back up
Access Ports	Two 2" Diameter
Interior Outlet	115V Duplex 3A
Doors	Two hinged glass, locking
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Six, locking
Shelves	Eight, Sliding/Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	2 x 1/5 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	12.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	525 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.

