

## SPECIAL OFFER # 62 - OVERSTOCK

### MODEL DH4-72GD-SS

## LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>72 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>2095 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>36°F to 46°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	75" x 33.75" x 81"
Interior [ Inches ]	72" x 27" x 55"
Exterior [ Centimeters ]	190 x 85 x 206
Interior [ Centimeters ]	183 x 69 x 140

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Stainless Steel
Interior	Stainless Steel
Insulation	Polyurethane Foam
Lockable Door	3 Doors, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	9 Adjustable Shelves
Shipping Weight	850 lbs.

### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   7.2A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

**Special Price Each: \$10,995.00**

**Price does not include freight costs.**

**Lead Time:** Approximately 1 weeks.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 3 Years Parts & Labor.  
 Additional 4 Years on Compressor Parts, No Labor.

**SPECIAL OFFER – SA62**