

MODEL DH29-74SD | SA66

LAB AND PHARMACY FREEZER -18°C TO -25°C

Capacity [Cu. Ft.]	74 Cu. Ft.
Capacity [Liters]	2095 Liters
Temperature Range [°C]	-18°C to -25°C
Temperature Range [°F]	0°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	78" x 34" x 81.5"
Interior [Inches]	74" x 28" x 60"
Exterior [Centimeters]	198 x 86 x 207
Interior [Centimeters]	188 x 71 x 152

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Timed Automatic
Refrigerant	R-404A
Exterior	Solid Steel
Interior	Metal
Insulation	Polyurethane Foam
Lockable Door	3 Doors, Keyed Lock
Door Material	Solid
Interior Light	2 Lights
Porthole	1/2" Diameter
Casters	4, Swivel
Storage	9 Adjustable Shelves
Shipping Weight	900 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	20 Amp Dedicated 8.4A
NEMA Configuration	NEMA 5-20

Controller

Display Type	Microprocessor, Digital
	Min / Max Memory
Alarms	High/Low Temp
Alarm Relay	Dry Contacts
	Power Failure Relay

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Years Parts & Labor.