

MODEL DHK4-10SD-K | SA72

LABORATORY GRADE REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	10 Cu. Ft.
Capacity [Liters]	283 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F

Dimensions	Width x Depth x Height
Exterior [Inches]	24" x 22" x 60"
Interior [Inches]	19.5" x 17" x 50"
Exterior [Centimeters]	60 x 55 x 152
Interior [Centimeters]	49 x 43 x 127

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle / Auto
Refrigerant	R600
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Leveling Legs	4, Adjustable
Storage	5 Adjustable Shelves
Shipping Weight	300 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 1.8A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor.