

MODEL LRT30W1HC | SA76

SOLID DOOR LAB REFRIGERATOR -2°C TO -16°C

Capacity [Cu. Ft.]	29.1 Cubic Feet
Capacity [Liters]	831 Liters
Temperature Range [°C]	-2°C to 16°C
Temperature Range [°F]	28°F to 61°F

Dimensions	Width x Depth x Height
Exterior [Inches]	31.25" x 38.75" x 85.25"
Interior [Inches]	24.5" x 31.75" x 59"
Exterior [Centimeters]	79 x 98 x 216
Interior [Centimeters]	62 x 80 x 150

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R290
Exterior	Steel, Painted White
Interior	Stainless Steel
Insulation	Polyurethane Foam
Interior Lighting	LED Ceiling Light
Self-Closing Door	1 Door, Automatic, Recessed
Lockable Door	Keyed Lock
Door Material	Insulated Steel
Porthole	1/2" Diameter
Casters	2 Swivel Locking, 2 Fixed
Shelving	4 Adjustable Shelves
Shelving Weight Limit	60 Lbs.
Approximate BTU/Hour	1465
Approximate Decibels	62 dB

Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 5A
NEMA Configuration	NEMA 5-15
Certification(s)	CSA Listed Energy Star Rated

Control:

Temperature Control	Microprocessor, Digital
Alarms	High / Low Temp
Min / Max Temperature	Display and Reset
Alarm Signal	Audible and Visual
Alarm Relay	Dry Contacts
Alarm Forwarding	Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor