

## MODEL DHH4-30SDFMS | SA80

# FLAMMABLE MATERIAL STORAGE REFRIGERATOR 1°C TO 10°C

<b>Capacity [ Cu. Ft. ]</b>	<b>30 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>849 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>1°C to 10°C</b>
<b>Temperature Range [ °F ]</b>	<b>33°F to 50°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	35.75" x 35.25" x 72.25"
Interior [ Inches ]	32" x 22.5" x 65"
Exterior [ Centimeters ]	90 x 90 x 184
Interior [ Centimeters ]	81 x 57 x 165

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R404A
Exterior	Painted Steel
Interior	Painted Steel
Insulation	High Density Urethane
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Porthole	3/4" Diameter
Leveling Legs	4
Storage	5 Adjustable Shelves
Shipping Weight	380 lbs.

### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   6A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

### Controller

Display Type	Mechanical Thermostat
--------------	-----------------------

### Digital Control Option\* \$ 1,250.00

\*Digital control option includes control & alarm system.

Note: If a temperature recorder is purchased, the digital control option must also be purchased.

» **Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.**

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 1 Year Parts & Labor.