

MODEL DHH4-14SDRX | SA81

EXPLOSION PROOF REFRIGERATOR 1°C TO 10°C

| | |
|---------------------------------|---------------------|
| Capacity [Cu. Ft.] | 14 Cu. Ft. |
| Capacity [Liters] | 396 Liters |
| Temperature Range [°C] | 1°C to 10°C |
| Temperature Range [°F] | 33°F to 50°F |

| | |
|--------------------------|------------------------|
| Dimensions | Width x Depth x Height |
| Exterior [Inches] | 23.5" x 24.5" x 73" |
| Interior [Inches] | 60 x 62 x 185 |
| Exterior [Centimeters] | 19.75" x 16.75" x 63" |
| Interior [Centimeters] | 50 x 42 x 160 |

| | |
|-----------------|-----------------------|
| Compressor | Hermetic Compressor |
| Air Circulation | Forced Air |
| Defrost Method | Cycle |
| Refrigerant | R600a |
| Exterior | Painted Steel |
| Interior | Painted Steel |
| Insulation | High Density Urethane |
| Lockable Door | 1 Door, Keyed Lock |
| Door Material | Insulated Steel |
| Leveling Legs | 4 |
| Storage | 5 Adjustable Shelves |
| Shipping Weight | 200 lbs. |

Electrical

| | |
|-------------------------|----------------------------|
| Voltage: | 115 Volts 60 Hertz 1 Phase |
| Amperage Line Running | 15 Amp Dedicated 3A |
| Electrical Plug | Not Supplied (Hardwired) |
| Certification(s) | ETL Listed |

Controller

| | |
|-----------------------|-------------------------|
| Display Type | Microprocessor, Digital |
| Alarms | High/Low Temp |
| Min / Max Temperature | Display and Reset |
| Alarm Relay | Dry Contacts |

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines for Hazardous Locations Class 1, Division II, Group C & D**

Warranty 2 Years Parts & Labor.