

SO-LOW

Environmental Equipment Co., Inc.

10310 Spartan Dr., Cincinnati, Ohio 45215 USA

(513) 772-9410 e-mail: sales@so-low.com

Fax (513) 772-0570 http://www.so-low.com

MODEL DHH4-30SDRX | SA82

EXPLOSION PROOF REFRIGERATOR 1°C TO 10°C

Capacity [Cu. Ft.]	30 Cu. Ft.
Capacity [Liters]	849 Liters
Temperature Range [°C]	1°C to 10°C
Temperature Range [°F]	33°F to 50°F

Dimensions	Width x Depth x Height
Exterior [Inches]	35.75" x 35" x 77.25"
Interior [Inches]	32" x 22.5" x 65"
Exterior [Centimeters]	90 x 88 x 196
Interior [Centimeters]	81 x 57 x 165

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R404A
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Leveling Legs	4
Storage	5 Adjustable Shelves
Shipping Weight	390 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6A
Electrical Plug	Not Supplied (Hardwired)
Certification(s)	UL Listed

Controller

Display Type	Mechanical Thermostat
--------------	-----------------------

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines for Hazardous Locations Class 1, Division II, Group C & D**

Warranty 1 Year Parts & Labor.

Specifications subject to change without notice